

UriCap Female Data Sheet

General Sales information

UriCap Female Devices and Extension Tubes can be ordered in bulk or packaged in boxes according to the product list below.



US pending patent application: 16/472086

FDA Registered Class I Medical Device

Components:

The UriCap female is assembled of two parts:

UriCap Female Device specifications









Name	UriCap Female device
PN	FPR02020L
Material	Medical grade silicone
Size packed	15*10*2.5 cm 6*4*1 in
Weight packed	7 gr 0.24 oz
Shelf life	3 years
Storage conditions	10-30°C/ 50-86°F
US Reimbursement code (HCPCS)	A4328
Comments	Single Use device/24hrs

Extension Tube specifications



Name	Extension tube
PN	FPR00001
Material	Tube- Medical grade Latex Connector- Medical grade polypropylene
Size packed	25*7.5*1.5 cm 9.8*3*0.6 in
Weight packed	8.3 gr 0.3 oz
Shelf life	3 years
Storage conditions	10-30°C/ 50-86°F
US Reimbursement code (HCPCS)	A4331
Comments	Contains Latex The Tube should be replaced with the urine bag

Product List:

	Part #	Description	Weight	Size
	FPR02020L	1 UriCap Female device	7 gr 0.24 oz	15*10*2.5 cm 6*4*1 in
	FPR00001	1 Extension Tube	8.3 gr 0.3 oz	25*7.5*1.5 cm 9.8*3*0.6 in
	KIT00005-10	10 UriCap Female sets in box (10 Devices + 10 Tubes)	230 gr 8.11 oz	24*14*9 cm 9.4*5.5*3.5 in
	FPR02020L-20	20 UriCap Female devices refill in box	202 gr 7.12 oz	24*14*9 cm 9.4*5.5*3.5 in
	FPR00001-20	20 Extension Tubes refill in box	227 gr 8 oz	24*14*9 cm 9.4*5.5*3.5 in
	KIT00010	UriCap Female Demonstration kit	334 gr 11.78 oz	22.5*9.7*5.5 cm 8.8*3.8*2.2 in

Local distributor

TillaCare Ltd.

2 Etgar st., 390321 Tirat
Carmel, Israel
E-mail: sales@tillacare.com



for more information

www.tillacare.com/training



TillaCare's mission is to bring to market cutting edge urine collection for incontinent men & women.

TillaCare's first product, UriCap Female, services a large need among the growing population of elderly women. We believe it will not only increase their dignity and health, but also decrease health care costs by preventing skin complications, reducing Urinary Tract Infections and the use of antibiotics.

UriCap is a an external, leak-free urine collection device designed to fit the female anatomy. It connects to a standard urine drainage bag.

UriCap is particularly useful for women that suffer from urinary incontinence with partial immobility, that are sitting up or are bedridden.

The device offers an innovative solution for patients in hospitals, nursing homes or home-care.



UriCap benefits

- Both leak-free and external
- Made of Medical Grade Silicone
- Takes a few moments to apply or remove
- No glue, tape, suction or other mechanical holders involved

For the patient

- An alternative to catheters and diapers
- Supports a full night sleep without diaper changes
- Comfortable and keeps the skin dry, day and night
- Prevents unpleasant odors and maintains hygiene
- Improves quality of life

For the caregiver

- Fast and easy to apply. No sterile procedure needed
- Exchanged only once daily.
- Less heavy lifting and improves work flow
- Enables non-invasive monitoring of urine
- Prevents excessive humidity, when treating burns, dermatitis and bed sores