

Tilla Care Product Catalogue

General Sales information

UriCap Female Devices and Extension Tubes SX-L can be ordered separately in bulk or packaged in kits and boxed according to the product list below.



US patent No.:
11,666,474

FDA Registered Class I
Medical Device

Components:

The UriCap Female is assembled from two parts, the UriCap Female and the Extension Tube.

Products in the catalogue and in the kits has been tested with the UriCap Female for optimum performance and functionality.

UriCap Female Device specifications

Extension Tube SX specifications



Name	UriCap Female Device	Name	Extension Tube SX-L
PN	FPR02020L	PN	FPR00004
Material	Medical grade silicone	Material	Tube- Medical grade silicone Connector- Medical grade polypropylene
Size packed	6*4*1 in. 15*10*2.5 cm.	Size packed	5*5.5*0.6 in. 25*7.5*1.5 cm.
Weight packed	0.24 oz 7g	Tube Length	11 in. 28 cm
Shelf life	3 years	Weight packed	0.33 oz 9.6 gr
Storage conditions	50-86°F 10-30°C	Shelf life	3 years
US Reimbursement code (HCPCS)	A4328	Storage conditions	50-86°F 10-30°C
Comments	Single Use device. certified for maximum 24 hrs. continuous use	US Reimbursement code (HCPCS)	A4331
		Comments	The Extension Tube should be replaced with the urine bag

Product List:

Image	Part #	Description	Quantities	Weight	Size
	FPR02020L	UriCap Female Device	Bulk is minimum order of 1000 pcs	247 oz 7 kg.	23*20.5*16 in. 58*52*41 cm.
	FPR00004	Extension Tube SX-L	Bulk is minimum order of 250 pcs	77 oz 2.5 kg.	19*14.5*11 in. 49*37*28 cm.

Local distributor:

Manufactured by:

TillaCare Ltd.


2 Etgar st., 390321 Tirat Carmel,
Israel

E-mail: sales@tillacare.com



For more information:

www.tillacare.com/training

Image	Part #	Description	Quantities	Weight	Size
	HRD00006	Securement Device	50	11.1 oz 315 gr.	7.3*7.3*3.7 in. 18.5*18.5*9.5 cm
	HRD00004	Lubricating Jelly	144	21.3 oz 604 gr.	3.3*5.5*5.1 in. 8.5*14*13 cm.
	KIT00010	UriCap Female Demonstration Kit	-	11.8 oz 335 gr.	8.8*3.8*2.2 in 24*14*9 cm
	HRD00011	Urine Bag Greystone	20	-	-

Kit Product List:

Image	Part #	Description	Kit Items	Quantities	Size
	KIT00013-30	Base Kit	1pc FPR02020L UriCap Female device 1pc FPR00004 Extension Tube SX-L 1pc HRD00004 Sterile Lubricating Jelly	30 kits 24.7 oz 700 g	10.5*8*7 in. 27*21*19 cm
	KIT00013-5H	Home Sample Kit (Box of 5)	5 KIT00013 Base Kits 2 HRD00006 Securement devices 1 HRD00008 Urine Bag Tilla 2L	5 kits 4.7 oz 135g	
	KIT00013-30H	Home Base Kit (30-Day supply)	30 KIT00013 Base Kits 10 HRD00006 Securement devices 2 HRD00008 Urine Bag Tilla 2L	30 kits 28.2 oz 800 g	10.5*8*7 in. 27*21*19 cm

	KIT00021	Information Pack	• DOC000050_B UriCap Female User Guide • IFU00003G UriCap Female SX IFU multilanguage • DOC000037 Paper Training Model		
--	----------	------------------	---	--	--